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Key words: cup method, desiccant method, water method, test dish, test atmosphere, seal

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The purpose of these tests is to obtain, by means of simple apparatus, reliable values of water vapor transfer through permeable and semipermeable materials, expressed in suitable units. These values are for use in design, manufacture, and marketing. (进行这些试验是为了通过简单的设备获得可信的可渗透材料以及半渗透材料的水蒸气透过值。这些测试值将用于设计、生产以及商贸。)

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ASTM E 96 将增重法和减重法都视为基础测试方法。The test methods are limited to specimens not over $1\frac{1}{4}$ in. (32mm) in thickness.....Two basic methods, the Desiccant Method and the Water Method, are provided for the measurement of permeance. (这两种试验方法可用于检测厚度不大于 $1\frac{1}{4}$ in. (32mm) 的试样……增重法和减重法两种基本方法都可以用来检测试样的水蒸气透过系数。)

In the Desiccant Method the test specimen is sealed to the open mouth of a test dish containing a desiccant, and the assembly placed in a controlled atmosphere. Periodic weighings determine the rate of water vapor movement through the specimen into the desiccant. (在增重法中, 试样被密封在装有干燥剂的透湿杯杯口处, 然后将整套透湿杯组件放在控制环境中。对透湿杯进行周期性称重来确定透过试样被干燥剂吸收的水蒸气透过率。) 增重法测试原理可参见图1。

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In the Water Method, the dish contains distilled water, and the weighings determine the rate of vapor movement through the specimen from the water to the controlled atmosphere. (在减重法中, 透湿杯中盛有蒸馏水, 通过称重数据确定从蒸馏水一侧渗透通过试样进入控制环境中的水蒸气透过率。) 图2 为减重法测试原理图。

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The vapor pressure difference is nominally the same in both methods except in the variation, with extremes of humidity on opposite sides. (两种测试方法的水蒸气压差相同, 都是试样一侧极湿而另一侧极干, 只是水蒸气压差相对于试样的方向有变化。)

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The test dish shall be of any noncorroding material, impermeable to water or water vapor. It may be of any shape. Light weight is desirable. A large, shallow dish is preferred, but its size and weight are limited when an analytical balance is chosen to detect small weight changes. (透湿杯应采用任何不腐蚀的、并且水或水蒸气无法渗透的材料。它可以是任意的形状。最好重量轻。首选透湿杯是大而浅的杯子, 但是选定了用于检测微小重量变化的分析天平后它的尺寸和重量也就受到了一定的限制。) **Different depths may be used for the Desiccant Method and Water Method, but a $\frac{3}{4}$ in. (19mm) depth (below the mouth) is satisfactory for**

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The material shall be sampled in accordance with standard methods of sampling applicable to the material under test. The sample shall be of uniform thickness. If the material is of nonsymmetrical construction, the two faces shall be designated by distinguishing marks. (应按照当前试样的相关取样标准进行材料试样的制取。试样应厚度均一。如果试样为非对称结构, 应使用指定符号对试样两面加以区分。) **Test specimen shall be representative of the material tested.** (测试试样应该能够代表当前测试的材料。) 在 ASTM E 96 中对在不同情况下的试样制取数目以及试验方式都有详细规定。

Attach the specimen to the dish by sealing (and clamping if desired) in such a manner that the dish mouth defines the area of the specimen exposed to the vapor pressure in the dish. Thoroughly seal the edges of the specimen to prevent the passage of vapor into, or out of, or around the specimen edges or any portion thereof. (通过密封 (如果需要可用夹紧方式) 将试样与透湿杯装夹在一起, 这样透湿杯的杯口就确定了试样的试验面积。完全密封试样的边缘以防止水蒸气进入、离开、在试样边缘周围或它的任何部分。)

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密封圈密封的适应性、可操作性要比蜡封好的多。**Gasketed types of seals are also in use on appropriately designed dishes. These simplify the mounting of the specimen, but must be used with caution, since the possibility of edge leakage is greater with gasketed seals than with wax seals.** (使用密封圈的密封方式在一些设计适当的透湿杯中也有应用。这简化了试样的装备, 但应小心使用, 因为使用密封圈的密封方式出现边缘泄漏的可能性比使用密封蜡要大。)

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Fill the test dish with desiccant.....Leave enough space so that shaking of the dish, which must be done at each weighing, will mix the desiccant. Attach the specimen to the dish and place it in the controlled chamber.....Weigh the dish assembly periodically, often enough to provide eight or ten data points during the test. At first the weight may change rapidly; later a steady state will be reached where the rate of change is substantially constant. (向透湿杯中添加干燥剂……留出足够的空间以满足摇动透湿杯混合杯内干燥剂的目的, 这一步是每次称重都必须要做的。将试样装夹在透湿杯上, 并将装夹有试样的透湿杯放进测试腔中……对透湿杯进行周期性称重, 需要在整个试验中获得 8 到 10 组数据。一开始重量可能会迅速变化; 之后会达到一个(渗透)稳定的状态, 这时的变化率就非常稳定了。)

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Fill the test dish with distilled water.....The air space thus allowed has a small vapor resistance, but it is necessary in order to reduce the risk of water touching the specimen when the dish is handled. Such contact invalidates a test on some materials such as paper, wood, or other hygroscopic materials. Attach the specimen to the dish. Weigh the dish assembly and place it in the controlled chamber on a true horizontal surface. (向透湿杯中添加蒸馏水……留出的空气隙可有很小的水蒸气阻抗性能, 但是空气隙的存在可减小在处理透湿杯时水接触到试样的风险。这种接触可以使得对一些材料的试验直接无效, 例如纸张、木材、以及其它吸湿性材料。将试样装夹在透湿杯上。称量试样、透湿杯组合的重量并将它放置在测试腔中绝对水平的表面上。) 之后的试验过程与增重法完全一致: 对透湿杯进行周期性称重, 直至重量变化率保持恒定。

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Calculate the water vapor transmission, and permeance as follow. (按照下式计算水蒸气透过率和水蒸气透过系数。)

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G/t = slope of the straight line, g/h (直线段的斜率, 单位是 g/h)

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When the specimen area is larger than the mouth area, this overlay upon the ledge is a source of error, particularly for thick specimens. The overlay material results in a positive error, indicating excessive water vapor transmission. The magnitude of the error is a complex function of the thickness, ledge width, mouth area, and possibly the permeability.

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6.4 dummy 试样

When results of water vapor transmission are expected to be less than 0.05 perm, a dummy specimen is strongly recommended. Such a dummy specimen should be attached to an empty cup in the normal manner. The environmental effects of temperature variation and buoyancy variability due to barometric pressure fluctuation can be arithmetically tared out of the weighing values. This precaution permits earlier and more reliable achievement of equilibrium conditions. (如果预计试样的水蒸气透过系数低于 **0.05 perm** (本标准中 **1 perm** 是 $5.7 \times 10^{-11} \text{ kg/Pa}\cdot\text{s}\cdot\text{m}^2$), 强烈推荐使用 **dummy** 试样。将试样按照通常的装夹方法装夹在空透湿杯上。温度变化以及由于大气压变化引起的浮力变化的这些环境影响可以在称重数据中去除掉。这样做可使得判断进入平衡状态的时间更早而且更加可靠。)**As a further precaution when gasketed seals are used instead of preferred sealants, a blank test run is suggested using glass or metal as a dummy specimen.** (当使用密封圈来代替首选的密封剂, 建议使用玻璃或是金属作为 **dummy** 试样进行空白试验。)

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